

DIP SWITCH SETTINGS (- = Irrelevant)	ON	OFF
DS/W1 in main pieceboard 501335	configuration	normal
DIP 1 - dental unit mode	disabled	enabled
DIP 2 - hot water outlet for spraying on CN16	-	-
DIP 3 - not used	-	-
DIP 4 - not used	-	-

ADJUSTMENTS

Main pieceboard 501335	
TR1	Optic fibre light intensity
TR2	Minimum micromotor speed
TR3	Maximum micromotor speed (Max voltage scaler)
TR4 e TR5	Free
Micromotor	pieceboard 50 rev. 501046 and 900 rev. 501262
P1	Analogue-to-digital switching
P2	Micromotor torque
Pushbutton	pieceboard 501412
TR1	Display contrast adjustment




WARNING SIGNS

Main pieceboard 501335	
DL1	Free
DL2	LED - enable delivery unit 4
DL3	LED - enable delivery unit 3
DL4	LED - enable delivery unit 2
DL5	LED - enable delivery unit 1
DL6	LED - enable optic fibre
DL7	LED - SBPO out
DL8	LED - enable tooth scaler out
DL9	LED - analog exchange MM or tooth scaler out
DL10	LED - inlet of delivery unit 4
DL11	LED - inlet of delivery unit 3
DL12	LED - inlet of delivery unit 2
DL13	LED - inlet of delivery unit 1
DL14	LED on during EEPROM re-initialization
DL15	LED - 28V's power on
DL16	LED - VCC 5V power on
DL17	LED - 24VAC2 power on
DL18	LED - enable optic fibre
DL19	LED in SBPO monitor
DL20	LED - Reset for prog.
DL21	LED - Press for prog.
DL22	LED in mot. safety dev.
DL23	LED - Invert MM out
DL24	LED - delivery unit interblock output




STORAGE DEVICE SETTINGS

In case of replacement of the general pcbboard, move U2 in the new pcbboard to keep the settings made.

MAIN PIECEBOARD 501335



DISPLAYED ERRORS ARE DESCRIBED IN THE INSTRUCTION HANDBOOK



REV. 4, UPDATED ON 17/03/2017

OMIS
CARLE & SANZANO
 (PUNTO) -DENT-
 OMS reserves all property rights on this document and does not allow reproduction or distribution to third parties without prior written consent.

DESCRIPTION	GENERAL DIAGRAM PUNTO KART NEW WITH UNIVERSAL TOP (TABLE 1/5)
MATERIAL	
DATE	30/05/2012
DRAWN BY	
UPPER SECT.	
MAT. CODE	
PRL. CODE	
SCALE	
DWG. NO.	520549

WARNING SIGNS

Assistant table pieceboard 501371	DL1	LED - press for flashing
	DL2	LED reset ext per cel.
	DL3	LED - ASP camillas 1
	DL4	LED - ASP camillas 2
	DL5	LED - enable camera out
	DL6	LED flashing with LIN bus active
	DL7	LED flashing if program is being correctly executed
Water unit pieceboard 501372	DL1	LED - press for flashing
	DL2	LED reset ext per cel
	DL3	LED - full level
	DL4	LED - intermediate level
	DL5	LED flashing with LIN bus active
	DL6	LED flashing if program is being correctly executed
	DL7	LED - enable camera
	DL8	LED - operating lamp out
	DL9	LED - spittoon out
	DL10	LED - cup 2 out
	DL11	LED - cup 1 out
	DL12	LED - CA/AP out
	DL13	LED - washing S.V. out
	DL14	LED - dental chair release S.V. out
	DL15	LED - part out
	DL16	LED - asp relay out
	DL17	LED - partialization S.V. out
	DL18	LED - water exchange S.V. out
	DL19	LED - part 2 out
Transformer pieceboard 501373	DL1	LED - fuse F4 in lamp malfunctioning
	DL2	LED - 28VAC fuse F5 in lamp malfunctioning
	DL3	LED - 24VAC2 fuse F6 in lamp malfunctioning
	DL4	LED - 28VS fuse F8 in lamp malfunctioning
	DL5	LED - 24VACT fuse F7 in lamp malfunctioning
	DL6	LED - VS fuse F9 in lamp malfunctioning
	DL7	LED - 28VDC power
	DL8	LED - CA/AP relay out monitor
	DL9	LED - ASP relay out monitor
	DL10	LED flashing if program is being correctly executed
	DL11	LED press for flashing
	DL12	LED reset ext per cel
	DL13	LED flashing with LIN bus active
Foot control pieceboard 501343	DL1	LED - inputs active
	DL2	LED - BUS active


JUMPER AND ADJUSTMENTS

Water unit pieceboard 501372		
Jp4, Jp5	1-2 closed: N.O. contact is obtained in E/C PAR2 (LIM.2) of CN6	
	2-3 closed: 28Vdc is obtained in E/C PAR2 (LIM.2) of CN6	
Jp6, Jp7	1-2 closed: N.O. contact is obtained in E/C PAR1 (LIM.1) of CN6	
	2-3 closed: 28Vdc is obtained in E/C PAR1 (LIM.1) of CN6	

DIP-SWITCH SETTINGS (- = irrelevant)

DSW1 in water unit pieceboard 501372		
DIP 1 - enable spittoon S.V. after cup filling	ON	enabled
DIP 2 - control of spittoon	on/off	timed
DIP 3 - delay in aspiration switch-off	0 s	5 s
DIP 4 - free	-	-

CUP AND SPITTOON TIME SETUP
 Water group pboard 501372
 Hold the button pressed until a beep is heard. Release the button when the cup has been filled to the desired level or the required time for basin rinsing has elapsed.



REV. 4 UPDATED ON 17/03/2017

OMIS
 OSILE IN SERVIZIO
 (90000) - 0431-24111

NOTE: All dimensions are in millimeters unless otherwise specified. All drawings are subject to change without notice. All dimensions are subject to change without notice. All dimensions are subject to change without notice.

DESCRIPTION	GENERAL DIAGRAM PUNTO KART NEW WITH UNIVERSAL TOP (TABLE 2/5)
MATERIAL	
DATE	30/05/2012
DRAWN BY	
UPPER SECT.	
MAT. CODE	
PRJ. CODE	
SCALE	
DWG. NO.	52054,9

FUSES



F	8A1 - general (230 Vac)
F1	Transformer pieceboard 501373
F2	4A1 - primary transformer fuse (230 Vac)
F3	5A1 - ASP OUT (230 Vac)
F4	5A1 - dent al chair motors (230 Vac)
F5	6.3A1 - operating lamp (16 or 22.8 Vac)
F6	8A1 - pboard and bus power supply 24VS (28 Vac)
F7	6.3A1 - syringes, 24 Vac motors, polym. lamp (24 Vac)
F8	6.3A1 - boiler, separator, ca S.V. (24 Vac)
F9	6.3A1 - prog. turbine, tooth sealer, cc S.V. (Vad)

WARNING SIGNS

DL1	LED - failure fuse F4 lamp
DL2	LED - failure fuse F5 28VAC
DL3	LED - failure fuse F6 24VAC2
DL4	LED - failure fuse F8 28VS
DL5	LED - failure fuse F7 24VAC1
DL6	LED - failure fuse F9 VS
DL7	LED - 28VDC power
DL8	LED - CA/AP relay out monitor
DL9	LED - ASP relay out monitor
DL10	LED flashing if program is being correctly executed
DL11	LED press for flashing
DL12	LED reset ext per col
DL13	LED flashing with LIN bus active
Foot control pieceboard 501343	
DL1	LED - inputs active
DL2	LED - BUS active

WARNING:
INSERT A FERRITE ON THE CABLES OUTGOING FROM CH3 AND CH5 NEAR THE PIECEBOARD MOD. "FERRITE K43194-8P" cod. 438449

WARNING:
INSERT A FERRITE ON THE CABLES OUTGOING FROM HS AND HT NEAR THE PIECEBOARD MOD. "FERRITE K43194-8P" cod. 438449



DIP-SWITCH SETPOINT-S:FOOT CONTROL FUNCTIONS

DSW in Foot control pieceboard 501343	1	2	3	4
Levers - controlled dental chair movements	on	off	off	off
Levers - controlled dental chair movements, main lever - controlled cups and operating lamp	on	on	off	off

REV. 4 UPDATED ON 17/03/2017

OMIS		GENERAL DIAGRAM PUNTO KART NEW WITH UNIVERSAL TOP (TABLE 3/5)	
DESCRIPTION	MATERIAL	UPPER SECT.	
DATE	30/05/2012	MAT. CODE	
DRAWN BY		PRJ. CODE	
SCALE		DWG. NO.	520549


OMIS reserves all property rights on this drawing and does not allow reproduction or partial reproduction without its authorization.

OMIS IN SERVIZIO (PRODOTTO) - ITALY -

LIST OF DISPLAYED ERRORS

E000	foot control pboard is not connected or is malfunctioning
E001	pushbutton panel pboard is not connected or is malfunctioning
E002	pboard of the floor-mounted box not connected or malfunctioning
E003	pboard of the assistant table is not connected or is malfunctioning
E004	dental chair pboard is not connected or is malfunctioning
E005	main pboard is not connected or is malfunctioning
E006	water group pboard is not connected or is malfunctioning
E007	faulty drive micromotor
E016	other error
E017	power protection
E018	motor phase missing
E019	serial communication timeout
E020	invalid EEPROM error
E021	over temperature error
E022	under voltage error
E023	over voltage error
E024	additional board disconnected (after startup) error
E025	timeout with additional board error
E026	invalid additional board version error
E027	invalid EEPROM in additional board error
E028	invalid mode error
E029	not used (free)
E030	frame error
E031	communication error between RS232
E032	28DC power fuse broken (fuse F8)
E033	24AC power fuse 1 broken (fuse F7)
E034	Lamp fuse broken (fuse F4)
E035	24AC power fuse 2 broken (fuse F6)
E036	output suction short circuit
E037	output assistant call short circuit
E038	overcurrent in branch 2, 24 AC
E039	pboards and bus overcurrent
E040	switching overheated
E041	stabilised 28 DC voltage error
E042	rectified 28 DC voltage error
E043	24 AC voltage error
E044	Lamp voltage error
E045	stabilised 24 DC voltage error
E046	overcurrent in branch 1, 24 AC
E047	lamp overcurrent
E064	"motor safety device" signal error
E065	EEPROM error water group pboard
E066	EEPROM error main pboard
E067	EEPROM error pushbutton panel pboard
E068	a touch is stuck in the instrument table pushbutton panel
E069	a touch is stuck in the assistant table pushbutton panel
E070	DAC error: dc micromotor and scaler
E071	overheating micromotor driver error
E072	battery error peda wireless
EP01	power solenoid valve short circuit
EP03	short circuit power solenoid valve to grounding system
EA01	water solenoid valve short circuit
EA03	short circuit water solenoid valve to grounding system
EC01	chip blower solenoid valve short circuit
EC03	short circuit chip blower solenoid valve to grounding system
ES01	spray solenoid valve short circuit
ES03	short circuit spray solenoid valve to grounding system
EL01	wash output short circuit
EL03	short circuit wash output to grounding system
EZ01	shutter output short circuit
EZ03	short circuit shutter output to grounding system
EX01	shutter output 1 short circuit
EX03	short circuit shutter output 1 to grounding system
EY01	shutter output 2 short circuit
EY03	short circuit shutter output 2 to grounding system
EH01	assistant call / open door output short circuit
EH03	short circuit assistant call/open door output to grounding system
ER01	suction output short circuit
ER03	short circuit suction output to grounding system
EO01	short circuit in distilled / standard water exchange output
EO03	short circuit to grounding system of distilled/strd. water exchange out
EB01	chair release piston output short circuit
EB03	short circuit to grounding system of chair release piston output

REV. 4 UPDATED ON 17/03/2017

 OMIS <small>OSALE E SALVIZIUM</small> <small>(PNOVA) -TAK-</small> <small>OMS delivers all necessary systems on</small> <small>time and cost and allows reproduction or</small> <small>modification of the drawings without</small> <small>any restriction.</small>		DESCRIPTION		GENERAL DIAGRAM PUNTO KART NEW WITH UNIVERSAL TOP (TABLE 5/5)	
		MATERIAL		DATE	30/05/2012
DRAWN BY				PRJ. CODE	
SCALE				DWG. NO.	520549